



Maximum LOAD CAPACITY: 300 lbs.

Limited Lifetime Warranty

A-Riser and Bryan Justice make every effort to assure that our products meet high quality and durability standards, and warrants to the original purchaser that this product is free from defects in materials and workmanship for the life of the original purchaser. This warranty does not apply to damage due directly or indirectly, to misuse, abuse, negligence or accidents, repairs or alterations outside our facilities, criminal activity, improper installation, improper use, normal wear and tear, or to lack of maintenance. We, and all associates and affiliates, shall in no event be liable for death, injuries to persons or property, or for incidental, contingent, special or consequential damages arising from the use of our product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation of exclusion may not apply to you. THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS.

To take advantage of this warranty, the product or part must be returned to us with transportation charges prepaid. Proof of purchase date and an explanation of the complaint must accompany the merchandise. If our inspection verifies the defect, we will either repair or replace the product at our election or we may elect to refund the purchase price if we cannot quickly provide you with a replacement. We will return repaired products at our expense, but if we determine there is no defect, or that the defect resulted from causes not within the scope of our warranty, then you must bear the cost of returning the product.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state.



A-Riser Lift

Owner's Manual & Safety Instructions









**Read this material before using this product.
Failure to do so can result in serious injury.**

SAVE THIS MANUAL.


Congratulations on the purchase of your A-Riser Lift!

This Manual contains Important and Critical Safety and Operation information. Read the ENTIRE manual prior to assembling or using the A-Riser.

All of the information is important.

WARNING SYMBOLS AND DEFINITIONS	
	SAFETY ALERT is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
	DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.
	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	CAUTION indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
	NOTICE is used to address practices not related to physical injury.
	SAFETY INSTRUCTIONS (or equivalent) signs indicate specific safety-related instructions or procedures.

Important Safety Information

	Failure to heed these warnings may result in personal injury:
---	--

1. Study, understand, and follow all instructions before operating this device.
2. Do not exceed the rated capacity. **The A-Riser Lift is rated at 300 lbs.**
3. Use only on hard level surfaces.
4. Inspect before every use; do not use if parts are loose or damaged. Do not use if even a minor crack, bend or other damage/deformation appears. Always check all

4. One side at a time, lift the Angle Iron cross braces, collapsing the side legs one to another.
5. Secure each side with a strap.
6. Remove the two pins on the wheel side top braces.
7. Lay the device down.
8. Remove the pins from the other two top braces.
9. Pull the legs together.
10. Fasten a strap around legs just above the feet.
11. You can now carry, drag or wheel the A-Riser to storage in a closet or other convenient location.

For Transport or Storage:

1. Simply wheel the device out to your car and place the top in the floorboard of one side, and the legs stick out the opposing window or rest on the seat back, depending on the size of your vehicle.
2. OR toss it in the trunk. Depending on the size of your vehicle, close or strap the lid.
3. OR if you have an SUV or Van, simply slide the A-Riser in for transport.
4. Store the A-Riser in a closet or other location. Place top support on floor with legs up.



4. Begin lifting by cranking the winch handle slowly.



WARNING Monitor for any discomfort, stress or pain. If so, stop immediately and lower the person back to the floor.

5. If the person is comfortable, continue cranking, lifting the person until their seat is just high enough to clear the seat plank or wheel chair. Do not raise higher than necessary.
6. Place the Seat plank (handles down) or roll your wheel chair under them and lock the wheels.
7. With your hand on the winch crank handle, apply a bit of upward cranking, and flip the winch thumb switch to the Down position. **DO NOT RELEASE THE CRANK!**



WARNING DO NOT remove your hand from the crank handle. Maintain pressure and control to prevent freefall.

8. Lower the person to the seat and lower more for slack in the cable.
9. Disconnect the hook and place it in the top braces out of the way.
10. Unfasten the snaps from the lifting belt, remove the spreader bar, unbuckle the lifting Strap and remove.
11. The fallen person is now up, in a sitting position or in a wheel chair. The A-Riser is also well suited as an additional standing support aid.

Collapsing and storing the A-Riser

1. Once the lift is accomplished, and the non-riser is risen and safe, it is time to collapse the A-Riser and put it away.
2. The collapse process is simply the reverse of the set-up.
3. Remove the 2 pins from the main diagonal braces.

bolts and nuts for proper tightness. Pivoting bolts are to be snug, but not overtight that would prevent pivoting.

5. This product is not a toy. Do not allow children to play with or near this device.
6. Use as intended only.
7. Do not use the A-Riser when you are tired or under the influence of alcohol or drugs. Read warning labels on prescriptions to determine if your judgement or reflexes are impaired while taking drugs. If there is any doubt, do not operate this device.
8. Use only on a dry, level, flat, firm surface free of grease and debris. The surface must be capable of supporting the device and load at the four leg points.
9. Be aware, the device folds, and thus provides “Pinch points” Do not stick your fingers (or other appendages) in where they might get pinched.
10. **When lowering a person**, your hand **MUST** remain on the winch crank handle. Your hand and pressure must not release. Keep your hand on the crank handle with adequate pressure to prevent a free fall .
11. NEVER lift a person higher than needed to get their seat to the bench or wheel chair.
12. NEVER attempt to lift a person who is injured.
13. Broken ribs are serious. **This device is for use of uninjured persons only.**
14. If you are unsure as to the injured/uninjured status of the individual to be lifted, DO NOT LIFT. Call for professional medical assistance.
15. The A-Riser cable must be centered over the lift strap on the person. Do not attempt to lift at an angle.
16. Always use the spreader bar.
17. Be sure all connections are secure, latched, buckled or fastened as required.

ASSEMBLY INSTRUCTIONS

1. Remove the packing materials from the device.
2. Place the A-Riser mainframe on the floor, casters down.
3. Undo the fastener securing straps.
4. Spread the legs apart.
5. The winch is loose and rotated for shipping. Rotate the winch to align the two holes, insert the pins and fasten.
6. Place the hole in the handle over the threaded crank stud on the winch and screw on the nut. Use the 3/4" wrench to tighten the nut securely.

OPERATION



Read the ENTIRE IMPORTANT SAFETY INFORMATION section at the beginning of this document including all text under subheadings therein before set up or use of this product.

Prepare the fallen person for the lift:



WARNING Check the fallen person carefully for injury. If injured, **DO NOT** use the device. Call for medical assistance.

1. If the person is not injured, roll them on their side, place the lifting strap flat on the floor (with the D rings down) at the shoulder blades, perpendicular to their body.
2. Roll the person on their back on top of the lifting strap.
3. Fasten the strap with the buckle comfortably and securely. There should be a D ring at each armpit or slightly above.
4. Make sure the fallen person is comfortable.
5. Connect the spreader bar to the D rings with the Snap Links.



WARNING Be sure the snap latches are fully engaged.

6. Lay the spreader bar on their chest.

Set up of the A-Riser mainframe:

1. Place the A-Riser mainframe on the floor, casters down, winch up.
2. Remove fastening straps.
3. Spread the legs, and at each of the four top sliding braces, align the holes, insert pins and latch.
4. Raise the unit to a standing upright position, with the rubber feet on the floor.
5. One side at a time, spread the feet and legs while lowering the side brace.



CAUTION Avoid pinch points as brace is lowered.

6. Lower the brace securely into slot. Repeat for other side.
7. Align the holes in the sliding diagonal braces at the rear legs. Insert pins and latch.
8. The A-Riser is now easily moved by the operator standing with the winch on the right inside the unit. Lift using the two main diagonal braces and place the unit centered over the fallen person's chest.
9. Position the A-Riser so that all four feet are on the floor equally and completely stable.

Connecting the A-Riser and Lifting:

1. Flip the winch thumb switch to the Lowering position. While pulling down on the hook, crank the winch to lower the hook to the chest of the person on the floor.
2. Flip the winch thumb switch back to the Lifting position.
3. Prepare the Seat Plank or your wheel chair close by and ready to position once the person is lifted.